

CONCENTRATED DISINFECTANT-DETERGENT PRODUCTS

SURFANIOS PREMIUM

- » Proven efficiency against hospital strains.
- » Aldehyde free, chlorine free.
- » Wide compatibility materials/floor recovering.
- » Optimisation of ecotoxicological information.



INDICATIONS

Cleaning and disinfection of floors, walls, medical equipment and non-invasive medical devices.

COMPOSITION

N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, didecylmethylammonium chloride, excipients.

PACKAGING - REF.

500 doses, 20 ml each 1917.129
 12 dosing bottles, 1L each 1917.092
 4 cans, 5L + 20 ml dosing pump 1917.036

INSTRUCTIONS FOR USE



PRODUCT TO BE DILUTED
 i.e. 20ml for 8L

Active from 15 min.

ACTIVE AGAINST	ACCORDING TO STANDARDS
BACTERIA	EN 1040, EN 1276, EN 13727, EN 13697, Active against MDRB according to EN 1276, EN 13697
MYCOBACTERIA	Mycobacterium terrae: EN 14348, EN 14563
YEASTS	EN 1275, EN 13624, EN 1650, EN 13697
MOULDS	EN 1275, EN 1650, EN 13697
VIRUSES	According to EN 14476: HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Herpesvirus, Vaccinia virus, Coronavirus

SURFANIOS / SURFANIOS CITRON

- » Aldehyde free.
- » Product available since 1990.
- » Combination of detergent and disinfectant action in a simultaneous operation.



INDICATIONS

Cleaning and disinfection of floors, walls and medical equipment.

COMPOSITION

N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, didecylmethylammonium chloride, excipients.

PACKAGING - REF.

	SURFANIOS	SURFANIOS CITRON
500 doses, 20ml each	350.129	347.129
12 dosing bottles, 1L each	350.092	347.092
4 cans, 5L + 20 ml dosing pump	350.036	347.036

INSTRUCTIONS FOR USE



PRODUCT TO BE DILUTED
 i.e. 20ml for 8L

Active from 15 min.

ACTIVE AGAINST	ACCORDING TO STANDARDS
BACTERIA	EN 1040, EN 13727, EN 1276, EN 13697 Active against MDRB
MYCOBACTERIA	Mycobacterium terrae: EN 14348, EN 14563
YEAST	EN 1275, EN 1650, EN 13624
MOULDS	EN 1275, EN 13624
VIRUSES	Active against HIV-1, BVDV (surrogate of HCV), PRV (surrogate of HBV), Influenza virus A (H ₁ N ₁)

